

# Eye Morphology KMPCLA\_EYE\_001

## Purpose

To detect abnormalities in eye morphology.

## Experimental Design

- **Minimum number of animals** : 7M + 7F
- **Age at test:** Week 59
- **Sex:** We do not expect the results of this test to show sexual dimorphism

## Procedure

1. Examine the anterior of both eyes (e.g. with slit lamp) and record any abnormalities
2. Test the iris/pupil light response
3. Image abnormal eyes as a minimum or all eyes if capacity permits
4. Dilate both eyes
5. Examine the anterior and posterior of both dilated eyes (e.g. with slit lamp and ophthalmoscope) and record any abnormalities
6. Image abnormal eyes as a minimum or all eyes if capacity permits

OCT:

1. Turn on the OCT and start the database
2. Anaesthetize mouse
3. Prepare mouse eyes with drops and place contact lens (focal length 10 mm) on the right eye
4. Enter mouse data in the "Create new patient file" area and switch to the "Acquisition" window
5. Move the OCT camera to the right position and activate measurement mode
6. Place mouse collaterally to the OCT camera on the right side of a platform that is fixed in front of the OCT lens
7. Search the contact lens in the live picture of the fundus image field and place the pupil of the mouse eye in the centre of the window
8. Move the OCT camera such that OCT lens and contact lens touch each other
9. Focus the fundus picture by slightly moving up/down or forward/backward
10. Save fundus images
11. Set the "Ref.Arm" ruler such that the section of the retina is placed in the centre of the blue rectangle
12. Set the mode of measurement on "vertical, horizontal line"
13. Move the blue horizontal line in the fundus image field to the optic nerve level
14. Save images of retinal sections
15. Move the OCT camera to the left position

16. Repeat measurement procedure for the left eye

#### Scheimpflug Imaging:

1. Turn on the Pentacam and start the patient data management
2. Apply one drop 0.5% Atropine to each mouse eye for pupil dilation
3. Enter mouse data in the "Patient" group box and switch to the Scan menu
4. Activate the "1 Picture" modus in the "Image Options" area
5. Move Pentacam to the right position
6. Hold the mouse on a platform such that the vertical LED 475 nm light slit is orientated in the center of the right eye ball
7. Guarantee optimal focus by using the fine adjustment software tool in the adjustment window
8. Start imaging manually by pressing the "Start Scan" button
9. Scheimpflug images are saved automatically
10. Move Pentacam to the left position
11. Repeat measurement procedure for the left eye

## Notes

- As a minimum, all abnormalities should be imaged.
  - Where capacity permits, all mice can be imaged
- Majority of parameters can be analysed using the standard approach for assessing categorical data. To increase power for analysis purposes, where an abnormality is detected in the left, right or both eyes, the data may be combined to generate one "abnormal" category.
- Data for both eyes is recorded under one parameter to distinguish phenotypes of incomplete penetrance in individuals and if an observation for one or both eyes cannot be made, this is recorded as 'no data'. The IMPC analysis pipeline does not take into account whether an abnormality is fully penetrant or not and the same weight is given for an abnormal observations in one or both eyes. In cases where it is not possible to confirm if an abnormality is present or not, the data is not included in the statistical analysis. The following logic is applied in determining whether to include the data in analysis:
  - If at least one of the eyes shows an abnormality in a particular parameter, the data for that specimen will be included in the statistical analysis even if the other eye is marked as "no data".
  - If the eyes are marked as "no data", or one eye is normal and the other eye is "no data" for a particular parameter the data for that specimen will not be included in the statistical analysis.

## Data QC

Image QC is typically performed during data collection to ensure high quality images are captured whilst eyes are dilated etc.

## Parameters and Metadata

# Optical Coherence Tomography Equipment Model KMPCLA\_EYE\_039\_001 | v1.2

procedureMetadata

Req. Analysis: true      Req. Upload: false      Is Annotated: false

Options: Envisu R2200, Spectralis, EnvisuTM R-Series SDOIS,

---

## Pupil Dilation KMPCLA\_EYE\_013\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Options: both eyes dilated, right eye dilated, no data for both eyes, no data left eye, right eye dilated, left eye dilated, no data left eye, no data right eye, left eye dilated, normal, no data right eye,

---

## Vitreous KMPCLA\_EYE\_083\_001 | v1.1

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Options: right eye abnormal, no data for both eyes, normal, no data left eye, right eye abnormal, no data left eye, no data right eye, left eye abnormal, left eye abnormal, both eyes abnormal, no data right eye,

---

## Left anterior chamber depth KMPCLA\_EYE\_067\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

## Topical Anesthetic KMPCLA\_EYE\_044\_001 | v1.1

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

**Options:** Oxybuprocain, Atropine sulphate, Phenylephrine hydrochloride, Atropine, Hydrochloride, No anesthesia, Mydriacyl,

---

## Retina (combined) KMPCLA\_EYE\_092\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

### Derivation:

retinaCombined('KMPCLA\_EYE\_020\_001', 'KMPCLA\_EYE\_021\_001', 'KMPCLA\_EYE\_022\_001')

---

**VIP of right eye** KMPCLA\_EYE\_078\_001 | v1.1

seriesMediaParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

---

**Images Ophthalmoscopy** KMPCLA\_EYE\_050\_001 | v1.1

seriesMediaParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

---

**Right posterior chamber depth** KMPCLA\_EYE\_065\_001 | v1.2

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Unit Measured: um

---

**Ophthalmoscope Equipment Model** KMPCLA\_EYE\_035\_001 | v1.2

procedureMetadata

Req. Analysis: true      Req. Upload: false      Is Annotated: false

**Options:** Omega 180 / 60D, Micron III, OMEGA 180 / Superfield NC, SL4 4AA, Sigma 150K, Omega 500 Unplugged, Genesis, Xenon Nova 175W light source + HOPKINS optic 1218AA /Nikon D5100 + 85 mm f/1.8 lens, Genesis-DF, Genesis-D,

---

## Retinal Blood Vessels Structure KMPCLA\_EYE\_025\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** no data right eye, left eye abnormal, both eyes abnormal, right eye abnormal, normal, left eye abnormal, no data left eye, no data left eye, right eye abnormal, no data for both eyes, no data right eye,

---

## Scheimpflug Equipment Manufacturer KMPCLA\_EYE\_041\_001 | v1.4

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** false

**Is Annotated:** false

**Options:** Oculus GmbH,

---

## VIP of right fundus KMPCLA\_EYE\_074\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Optical Coherence Tomography Equipment ID KMPCLA\_EYE\_0

37\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Bulging eye KMPCLA\_EYE\_002\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data right eye, no data for both eyes, present both eyes,  
no data left eye, present right eye, no data left eye, absent, no data right eye, present left eye,  
present left eye, present right eye,

---

## Scheimpflug description KMPCLA\_EYE\_053\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

# Eyelid closure KMPCLA\_EYE\_005\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

**Options:** no data right eye, no data for both eyes, normal, no data left eye, right eye closed, no data right eye, left eye closed, both eyes closed, left eye closed, no data left eye, right eye closed,

-----

# Corneal opacity KMPCLA\_EYE\_008\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: true

**Options:** no data left eye, present left eye, no data left eye, present right eye, present right eye, no data right eye, no data for both eyes, no data right eye, present left eye, present both eyes, absent,

-----

# Corneal mineralization KMPCLA\_EYE\_084\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

**Options:** absent, present left eye, no data for both eyes, no data right eye, present left eye, no data left eye, present right eye, no data left eye, present right eye, present both eyes, no data right eye,

-----

# Fusion between cornea and lens KMPCLA\_EYE\_018\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data left eye, present both eyes, present right eye, no data right eye, no data left eye, present right eye, present left eye, no data right eye, present left eye, no data for both eyes, absent,

---

# Right inner nuclear layer KMPCLA\_EYE\_063\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

# Ophthalmoscope Lens Model KMPCLA\_EYE\_089\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

# Left vitreous humour thickness KMPCLA\_EYE\_088\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

-----

## Slit Lamp Equipment Model KMPCLA\_EYE\_032\_001 | v1.2

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

**Options:** SL 990, SL30, SL-15, SL 139, BQ 900 LED/IM-900, Micron III slit lamp extension, S350, SL-7E, SL130, 30 SL-M,

-----

## Scheimpflug Equipment ID KMPCLA\_EYE\_040\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

-----

## Retinal Blood Vessels KMPCLA\_EYE\_024\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

**Options:** no data right eye, no data right eye, left eye abnormal, no data left eye, right eye abnormal, no data left eye, right eye abnormal, no data for both eyes, left eye abnormal, both eyes abnormal, normal,

---

## Iris/Pupil KMPCLA\_EYE\_010\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data left eye, right eye abnormal, no data right eye, both eyes abnormal,  
no data for both eyes, right eye abnormal, no data left eye, normal,  
no data right eye, left eye abnormal, left eye abnormal,

---

## Right vitreous humor thickness KMPCLA\_EYE\_087\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

## Ophthalmoscope Equipment ID KMPCLA\_EYE\_033\_001 | v1.2

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

# Ophthalmoscope Equipment Manufacturer KMPCLA\_EYE\_034\_0

01 | v1.2

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

**Options:** Phoenix Research Labs, Heine / Volk, Phoenix, Karl Storz / Nikon, Keeler LTD, Heine, Haag-Streit, Kowa,

---

## Right anterior chamber depth KMPCLA\_EYE\_061\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

## Lens Opacity KMPCLA\_EYE\_017\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

**Options:** present both eyes, no data right eye, present left eye, absent, no data right eye, no data for both eyes, present left eye, no data left eye, present right eye, no data left eye, present right eye,

---

# Pupil Position KMPCLA\_EYE\_011\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

**Options:** right eye abnormal, both eyes abnormal, left eye abnormal,  
no data left eye, right eye abnormal, normal, no data left eye,  
no data right eye, left eye abnormal, no data right eye, no data for both eyes,

---

# Right corneal thickness KMPCLA\_EYE\_060\_001 | v1.2

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Unit Measured: um

---

# Left posterior chamber depth KMPCLA\_EYE\_071\_001 | v1.2

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Unit Measured: um

---

# Retina KMPCLA\_EYE\_020\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Options:** both eyes abnormal, no data left eye, right eye abnormal, left eye abnormal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, normal, no data right eye, no data for both eyes,

---

## Right total retinal thickness KMPCLA\_EYE\_062\_001 | v1.2

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Unit Measured:** um

---

## Mean left eye lens density KMPCLA\_EYE\_056\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Unit Measured:** %

---

## Max right eye lens density KMPCLA\_EYE\_058\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

Unit Measured: %

---

**Corneal ulcer** KMPCLA\_EYE\_085\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

**Options:** no data for both eyes, no data right eye, present both eyes,  
no data right eye, present left eye, no data left eye, present right eye, no data left eye,  
present left eye, absent, present right eye,

---

**Right eye diameter** KMPCLA\_EYE\_090\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Unit Measured: mm

---

**B-scan of left cornea and lens** KMPCLA\_EYE\_077\_001 | v1.1

seriesMediaParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

---

## Iris Pigmentation KMPCLA\_EYE\_015\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data for both eyes, both eyes abnormal, normal, no data right eye, left eye abnormal, no data left eye, right eye abnormal, no data left eye, right eye abnormal, no data right eye, left eye abnormal,

---

## Date Scheimpflug equipment last calibrated KMPCLA\_EYE\_048\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

## Eye Hemorrhage or Blood Presence KMPCLA\_EYE\_003\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data for both eyes, no data right eye, no data left eye, present right eye, absent, no data left eye, present right eye, no data right eye, present left eye, present both eyes, present left eye,

---

# Optical Coherence Tomography Equipment Manufacturer

KMPCLA\_EYE\_038\_001 | v1.2

procedureMetadata

Req. Analysis: true                      Req. Upload: false                      Is Annotated: false

Options: Heidelberg Engineering, Bioptigen,

---

## Left total retinal thickness KMPCLA\_EYE\_068\_001 | v1.2

simpleParameter

Req. Analysis: false                      Req. Upload: false                      Is Annotated: true

Unit Measured: um

---

## Mean right eye lens density KMPCLA\_EYE\_059\_001 | v1.1

simpleParameter

Req. Analysis: false                      Req. Upload: false                      Is Annotated: true

Unit Measured: %

---

## Date OCT equipment last calibrated KMPCLA\_EYE\_049\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Retinal Blood Vessels Pattern KMPCLA\_EYE\_026\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data left eye, normal, right eye abnormal, left eye abnormal,  
no data right eye, left eye abnormal, both eyes abnormal,  
no data left eye, right eye abnormal, no data right eye, no data for both eyes,

---

## VIP of left eye KMPCLA\_EYE\_079\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Slit Lamp observation KMPCLA\_EYE\_028\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Images Slit Lamp KMPCLA\_EYE\_051\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

**Slit Lamp Equipment ID** KMPCLA\_EYE\_030\_001 | v1.2

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

**Eyelid morphology** KMPCLA\_EYE\_004\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data right eye, right eye abnormal, no data left eye, normal,  
no data for both eyes, both eyes abnormal, no data right eye, left eye abnormal,  
left eye abnormal, no data left eye, right eye abnormal,

**B-scan of right retina** KMPCLA\_EYE\_072\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

# Slit Lamp Equipment Manufacturer KMPCLA\_EYE\_031\_001 | v1.2

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

Options: Kowa, CSO, Topcon, Haag-Streit, Phoenix Research Labs, MuLe, Zeiss,

---

## Eye KMPCLA\_EYE\_001\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: absent both eyes, present, absent left eye, absent right eye,

---

## Left corneal thickness KMPCLA\_EYE\_066\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

## Right outer nuclear layer KMPCLA\_EYE\_064\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

**Persistence of hyaloid vascular system** KMPCLA\_EYE\_027\_001 |

v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** present right eye, no data right eye, no data left eye,  
no data right eye, present left eye, absent, no data left eye, present right eye,  
no data for both eyes, present left eye, present both eyes,

---

**Retinal Pigmentation** KMPCLA\_EYE\_021\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

**Options:** no data right eye, left eye abnormal, normal, no data for both eyes,  
no data left eye, right eye abnormal, left eye abnormal, no data right eye, right eye abnormal,  
both eyes abnormal, no data left eye,

---

**Dilation Method** KMPCLA\_EYE\_043\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Options:** Atropine, None, Tropicamide+Phenylephrin, Atropine sulphate, Cyclopentolate hydrochloride, Tropicamide, Phenylephrine hydrochloride, Cyclopentolate hydrochloride+Phenylephrine hydrochloride,

---

## Ophthalmoscope Observation KMPCLA\_EYE\_029\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

## Min left eye lens density KMPCLA\_EYE\_054\_001 | v1.2

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Unit Measured:** %

---

## Scheimpflug Equipment Model KMPCLA\_EYE\_042\_001 | v1.4

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** false

**Is Annotated:** false

**Options:** Pentacam,

---

## Date Slit Lamp equipment last calibrated KMPCLA\_EYE\_046\_001

| v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Max left eye lens density KMPCLA\_EYE\_055\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

---

## Pupil Light Response KMPCLA\_EYE\_014\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data right eye, left eye abnormal, no data left eye, right eye abnormal, normal, left eye abnormal, no data left eye, right eye abnormal, no data right eye, no data for both eyes, both eyes abnormal,

---

## Corneal deposits KMPCLA\_EYE\_081\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data right eye, no data for both eyes, no data left eye, present right eye, present left eye, no data left eye, present both eyes, no data right eye, present left eye, present right eye, absent,

---

**B-scan of left retina** KMPCLA\_EYE\_073\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

**Synechia** KMPCLA\_EYE\_019\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** present both eyes, no data left eye, no data right eye, no data right eye, present left eye, no data left eye, present right eye, no data for both eyes, absent, present right eye, present left eye,

---

**Left inner nuclear layer** KMPCLA\_EYE\_069\_001 | v1.2

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Unit Measured:** um

---

## Experimenter ID KMPCLA\_EYE\_036\_001 | v1.1

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

---

## Corneal Sclerization KMPCLA\_EYE\_080\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** present both eyes, present right eye, no data left eye, present left eye, no data right eye, present left eye, no data right eye, no data for both eyes, absent, no data left eye, present right eye,

---

## General Anesthetic KMPCLA\_EYE\_045\_001 | v1.1

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Options:** Euthatal, No anesthesia, Ketamine+Medetomidine, Ketamine+Xylazine, Isoflurane, Avertin,

---

## VIP of left fundus KMPCLA\_EYE\_075\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Left outer nuclear layer KMPCLA\_EYE\_070\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

## Lens KMPCLA\_EYE\_016\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

**Options:** right eye abnormal, no data for both eyes, no data left eye, no data right eye,  
no data left eye, right eye abnormal, no data right eye, left eye abnormal,  
both eyes abnormal, left eye abnormal, normal,

---

## Pupil Shape KMPCLA\_EYE\_012\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** right eye abnormal, no data right eye, left eye abnormal, normal, both eyes abnormal, no data left eye, right eye abnormal, no data right eye, no data for both eyes, left eye abnormal, no data left eye,

---

## Min right eye lens density KMPCLA\_EYE\_057\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

---

## Date Ophthalmoscope equipment last calibrated KMPCLA\_EYE\_047\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Narrow eye opening KMPCLA\_EYE\_006\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** right eye abnormal, no data left eye, right eye abnormal, no data right eye, normal, no data left eye, no data right eye, left eye abnormal, both eyes abnormal, no data for both eyes, left eye abnormal,

---

**Retinal Structure** KMPCLA\_EYE\_022\_001 | v1.1

simpleParameter

**Req. Analysis:** false                      **Req. Upload:** true                      **Is Annotated:** false

**Options:** normal, right eye abnormal, no data right eye, left eye abnormal, no data right eye, no data for both eyes, both eyes abnormal, no data left eye, right eye abnormal, left eye abnormal, no data left eye,

---

**Optic Disc** KMPCLA\_EYE\_023\_001 | v1.0

simpleParameter

**Req. Analysis:** false                      **Req. Upload:** true                      **Is Annotated:** true

**Options:** no data right eye, left eye abnormal, right eye abnormal, no data left eye, both eyes abnormal, no data left eye, right eye abnormal, normal, no data for both eyes, left eye abnormal, no data right eye,

---

**Lacrimation** KMPCLA\_EYE\_086\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** absent, present left eye, present both eyes, present right eye,  
no data for both eyes, no data right eye, no data right eye, present left eye, no data left eye,  
no data left eye, present right eye,

---

**B-scan of right cornea and lens** KMPCLA\_EYE\_076\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

**Corneal vascularization** KMPCLA\_EYE\_009\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** absent, no data for both eyes, present left eye, present both eyes,  
present right eye, no data right eye, no data left eye, present right eye,  
no data right eye, present left eye, no data left eye,

---

**Iris transillumination** KMPCLA\_EYE\_082\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data left eye, right eye abnormal, no data right eye, left eye abnormal, both eyes abnormal, no data right eye, normal, right eye abnormal, no data left eye, no data for both eyes, left eye abnormal,

---

**Left eye diameter** KMPCLA\_EYE\_091\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** true

**Unit Measured:** mm

---

**Sheimpflug Lens description** KMPCLA\_EYE\_052\_001 | v1.1

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** false

---

**Cornea** KMPCLA\_EYE\_007\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** true

**Options:** no data right eye, left eye abnormal, right eye abnormal, no data left eye, right eye abnormal, no data for both eyes, normal, both eyes abnormal, no data left eye, no data right eye, left eye abnormal,

---

