

# X-ray JAXLA\_XRY\_001

## Purpose

Construct and analyse digital X-ray images in immobilised mice using a Faxitron X-Ray system or NTB digital X-ray scanner.

## Experimental Design

- **Minimum number of animals** : 4M + 4F
- **Age at test**: Week 75

## Procedure

Perform steps 3.1 and 3.2 unless the mouse is coming from DEXA analysis

3.1 Anesthetize the mice.

3.2 Monitor the animal carefully until unconsciousness by ensuring the mouse is adequately sedated.

3.3 Mount the anesthetized mouse onto an X-ray permeable plate. Curve the tail along the left side of the mouse. (Note: this will enable identification of the X-ray as a dorsal or ventral image.)

Faxitron X-ray system

3.3.1 Insert the plate into the Faxitron X-ray system.

3.3.2 Take an image using the software that is supplied by the X-ray scanner.

3.3.3 Set the machine to a voltage of 25kV with an integration time of 40ms.

3.3.4 Take an image using the iX-Pect software that is supplied by the X-ray scanner.

NTB digital X-ray scanner

3.4 Remove the mouse once the image is captured. Place the mouse on a heated mat, set at 37°C, in a cage and monitor closely until consciousness is fully regained.

3.5 Analyze the image.

3.6 Or continue with ABR.

# Notes

The protocol described herein employs anaesthesia using intraperitoneal (IP) administration; however, anaesthesia using isoflurane is also suitable.

[Digit integrity IMPC\\_XRY\\_031\\_001](#) *Definition:* The state of being whole and undivided.

## Parameters and Metadata

### Scanner equipment model JAXLA\_XRY\_046\_001 | v1.1

procedureMetadata

Req. Analysis: false      Req. Upload: false      Is Annotated: false

-----

### Number of thoracic vertebrae JAXLA\_XRY\_014\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

-----

### Free Ontology JAXLA\_XRY\_054\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

-----

### Humerus JAXLA\_XRY\_021\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: abnormal, imageOnly, normal, unobservable,

---

# XRay Images Skull Dorso Ventral Orientation JAXLA\_XRY\_051\_001 | v1.2

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

---

# Maxilla/Pre-maxilla JAXLA\_XRY\_003\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: abnormal, normal, unobservable, imageOnly,

---

# Joints JAXLA\_XRY\_027\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

**Options:** imageOnly, abnormal, unobservable, normal,

---

**Skull shape** JAXLA\_XRY\_001\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** true

**Options:** normal, imageOnly, abnormal, unobservable,

---

**Processes on vertebrae** JAXLA\_XRY\_020\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** true

**Options:** imageOnly, unobservable, normal, abnormal,

---

**Number of digits** JAXLA\_XRY\_028\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** true

---

## Transitional vertebrae JAXLA\_XRY\_060\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: imageOnly, yes, unobservable, no,

---

## Shape of ribcage JAXLA\_XRY\_059\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: imageOnly, normal, unobservable, abnormal,

---

## Date Xray equipment last calibrated JAXLA\_XRY\_047\_001 | v1.3

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Polysyndactylism JAXLA\_XRY\_062\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: imageOnly, unobservable, yes, no,

---

## Caudal vertebrae morphology JAXLA\_XRY\_074\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: normal, imageOnly, abnormal, unobservable,

---

## XRay Images Whole Body Lateral Orientation JAXLA\_XRY\_048\_001 | v1.2

seriesMediaParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Increments: Minimum 1

---

## Scapulae JAXLA\_XRY\_006\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: imageOnly, normal, unobservable, abnormal,

---

## Sacral processes JAXLA\_XRY\_066\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Options: unobservable, imageOnly, abnormal, normal,

---

## Number of ribs right JAXLA\_XRY\_008\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

---

## Fibula JAXLA\_XRY\_026\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: true

Options: normal, imageOnly, unobservable, abnormal,

---

## Scanner equipment ID JAXLA\_XRY\_044\_001 | v1.1

procedureMetadata

Req. Analysis: false      Req. Upload: false      Is Annotated: false

---

## Date Scanner equipment last calibrated

JAXLA\_XRY\_053\_001 | v1

.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Femur

JAXLA\_XRY\_024\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: imageOnly, unobservable, abnormal, normal,

---

## Fusion of ribs

JAXLA\_XRY\_011\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: imageOnly, unobservable, fused, not fused,

---

## Pelvis

JAXLA\_XRY\_012\_001 | v1.0

simpleParameter



**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** abnormal, unobservable, imageOnly, normal,

---

## **Kyphosis** JAXLA\_XRY\_057\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** no, imageOnly, yes, unobservable,

---

## **Number of cervical vertebrae** JAXLA\_XRY\_013\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

---

## **Tibia length** JAXLA\_XRY\_033\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Unit Measured:** mm

---

# Scanner equipment manufacturer JAXLA\_XRY\_045\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Tibia JAXLA\_XRY\_025\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: normal, imageOnly, abnormal, unobservable,

---

## Comment on XRay image JAXLA\_XRY\_032\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## XRay Images Hind Leg and Hip JAXLA\_XRY\_052\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

---

## Number of ribs left JAXLA\_XRY\_009\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

---

## Teeth JAXLA\_XRY\_005\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: abnormal, imageOnly, normal, unobservable,

---

## Mandibles JAXLA\_XRY\_004\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: abnormal, imageOnly, unobservable, normal,

---

## Fusion of processes JAXLA\_XRY\_061\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: unobservable, no, yes, imageOnly,

---

## XRay Images Whole Body Dorso Ventral JAXLA\_XRY\_034\_001 | v1.3

seriesMediaParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Increments: Minimum 1

---

## XRay Images Forepaw JAXLA\_XRY\_049\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

---

## Thoracic processes JAXLA\_XRY\_064\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** abnormal, imageOnly, normal, unobservable,

---

**Lumbar processes** JAXLA\_XRY\_065\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** unobservable, normal, imageOnly, abnormal,

---

**Ulna** JAXLA\_XRY\_023\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** abnormal, unobservable, normal, imageOnly,

---

**XRay Images Skull Lateral Orientation** JAXLA\_XRY\_050\_001 | v1.2

seriesMediaParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

**Increments:** Minimum 1

---

**Missing cranial rib** JAXLA\_XRY\_068\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

Options: Yes, No,

---

**Voltage settings** JAXLA\_XRY\_043\_001 | v1.1

procedureMetadata

Req. Analysis: false      Req. Upload: true      Is Annotated: false

---

**Number of lumbar vertebrae** JAXLA\_XRY\_015\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

---

**Radius** JAXLA\_XRY\_022\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: true

Options: normal, unobservable, abnormal, imageOnly,

---

## Shape of vertebrae JAXLA\_XRY\_018\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** unobservable, normal, imageOnly, abnormal,

---

## Shape of ribs JAXLA\_XRY\_010\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** abnormal, normal, unobservable, imageOnly,

---

## Shape of spine JAXLA\_XRY\_055\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** unobservable, imageOnly, normal, abnormal,

---

# Settings time of exposure JAXLA\_XRY\_039\_001 | v1.1

procedureMetadata

Req. Analysis: false      Req. Upload: true      Is Annotated: false

Unit Measured: s

---

# Cervical processes JAXLA\_XRY\_063\_001 | v1.1

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Options: unobservable, normal, imageOnly, abnormal,

---

# Settings level JAXLA\_XRY\_040\_001 | v1.1

procedureMetadata

Req. Analysis: false      Req. Upload: true      Is Annotated: false

---

# Pelvic vertebrae morphology JAXLA\_XRY\_073\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Options: unobservable, imageOnly, abnormal, normal,



---

## Lordosis JAXLA\_XRY\_058\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** yes, no, unobservable, imageOnly,

---

## Caudal processes JAXLA\_XRY\_067\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** abnormal, unobservable, imageOnly, normal,

---

## Experimenter ID (scan) JAXLA\_XRY\_041\_001 | v1.1

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

---

## Digit integrity JAXLA\_XRY\_031\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** yes, no, unobservable, imageOnly,

---

## **Zygomatic bone** JAXLA\_XRY\_002\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** abnormal, imageOnly, normal, unobservable,

---

## **Thoracic vertebrae morphology** JAXLA\_XRY\_071\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** normal, abnormal, unobservable, imageOnly,

---

## **Xray equipment model** JAXLA\_XRY\_038\_001 | v1.4

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** false

**Is Annotated:** false

**Options:** UltraFocus 100, MS-FX Pro, R\_mCT (Micro-CT), UltraFocus DXA, MX-20-DC4, MX-20, MX-20-DC12,

---

**Brachydactyly** JAXLA\_XRY\_030\_001 | v1.1

simpleParameter

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** true

**Options:** no, imageOnly, unobservable, yes,

---

**Clavicle** JAXLA\_XRY\_007\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** true

**Options:** abnormal, unobservable, normal, imageOnly,

---

**Number of caudal vertebrae** JAXLA\_XRY\_017\_001 | v1.1

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** true

---

## Rib morphology JAXLA\_XRY\_069\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: unobservable, normal, abnormal, imageOnly,

---

## Xray equipment manufacturer JAXLA\_XRY\_037\_001 | v1.4

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

Options: Faxitron Bioptics LLC, Bruker, Faxitron X-Ray Corporation,  
Rigaku Corporation (Micro-CT),

---

## Fusion of vertebrae JAXLA\_XRY\_019\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: no, imageOnly, unobservable, yes,

---

## Number of pelvic vertebrae JAXLA\_XRY\_016\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

---

## Syndactylism JAXLA\_XRY\_029\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: unobservable, imageOnly, Syndactylism, no syndactylism,

---

## Lumbar vertebrae morphology JAXLA\_XRY\_072\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: normal, unobservable, imageOnly, abnormal,

---

## Experimenter ID (analysis) JAXLA\_XRY\_042\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Alive JAXLA\_XRY\_035\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Cervical vertebrae morphology JAXLA\_XRY\_070\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: abnormal, unobservable, normal, imageOnly,

---

## Scoliosis JAXLA\_XRY\_056\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: yes, imageOnly, unobservable, no,

---

## Xray equipment ID JAXLA\_XRY\_036\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

