# Combined SHIRPA and Dysmorphology TC PIP\_CSD\_003

#### **Purpose**

SHIRPA and dysmorphology were originally always separate assessments. However they have recently been combined as assessments, so that they take place at the same time.

The purpose of the assessments is to examine mice for obvious physical characteristics, behaviors and morphological abnormalities.

Descriptions include abnormal locomotion/appearance/behavior/reflex reactions.

### **Experimental Design**

• Minimum number of animals: 7M + 7F

• Age at test: Week 26

• Sex: We would expect the results of this test to show sexual dimorphism

#### **Equipment**

- Viewing Jar
- SHIRPA arena
- Grid above arena
- Click Box
- Geotaxis grid
- Tube for contact righting
- Photo and/or video camera

#### **Procedure**

- 1. Allow the mice to acclimatise to the phenotyping room for a period of 30 minutes prior to testing.
- 2. Throughout the test note any vocalisation, aggression, salivation or unexpected behaviours.
- 3. During this test, note any unexpected dysmorphological characteristics (morphological irregularities are recorded by video sequences or photos).
- 4. Place the mouse in a clear cylinder over a wire grid and observe for activity and tremors.
- 5. Transfer the mouse out of the cylinder by removing the metal plate/grid whilst positioning 30cm over an arena and record the transfer arousal.

- 6. Record the number of  $10\text{cm}^2$  squares the mouse moves into in the first 30 seconds in the arena (locomotor activity).
- 7. Allow the mouse to move freely around the arena whilst being observed for gait and tail elevation.
- 8. Hold the click box approximately 30cm above the arena and press the button, record the response of the mouse.
- 9. Pick up the mouse by the tail and observe for limp grasping and trunk curl. Trunk curl must only be recorded if the mouse curls forward without twisting its body, bending to one side is not scored as a trunk curl.
- 10. Place the mouse in a small transparent tube. Turn the tube quickly so the mouse is fully upside down and record if the mouse rights itself.
- 11. Record any vocalisation and/or aggression which were observed throughout the entire test

#### **Notes**

- 1. If wiping down with ethanol prior to the use of equipment, make sure no ethanol residue remains as the ethanol may affect the behaviour of the animals.
- 2. The validity of results obtained from behavioural phenotyping is largely dependent on methods of animal husbandry. It is important that individuals following this procedure are experienced and aware of the animal's welfare, and is familiar with the animal being tested, in order to reduce the anxiety levels of the animal prior to testing.
- 3. The majority of mouse behavioural studies are age/sex/strain dependent. It is important to keep these parameters comparable throughout a single experiment.
- 4. Environmental factors may contribute to the levels of anxiety within the mouse. The temperature, humidity, ventilation, noise intensity and light intensity must be maintained at levels appropriate for mice. It is essential that the mice be kept in a uniform environment before and after testing to avoid anomalous results being obtained.
- 5. It is recommended that all phenotyping experimentation is conducted at approximately the same time of day because physiological and biochemical parameters change throughout the day.
- 6. When a number of mice are tested continuously, residual odours from the equipment used in the preceding test may affect the test results. The floor and walls of the arena, ruler, and metal net should be wiped clean before introducing the next mouse. To prevent infection, the equipment should be washed with water at the completion of the day's tests. Some specific pathogen-free facilities use ultraviolet irradiation when tests are not being performed. Care needs to be taken, however, to ensure that ultraviolet irradiation does not crack any acrylate equipment covered with residual alcohol.

#### **Parameters and Metadata**

Coat - color - head TCPIP\_CSD\_006\_001 | v1.0

simpleParameter

Options: No data, As expected, Not as expected,		
Coat - bair toyturo	Lannoaranco - hoad	TODID COD OUT OOM LINE O
simpleParameter	/ appearance - nead	TCPIP_CSD_015_001   v1.0
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: As expected, Not as	s expected, No data,	
Unexpected behave simpleParameter	viors TCPIP_CSD_037_001	v1.0
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: None, Jumping, Retropulsion, Other, Circling,		
Gait TCPIP_CSD_033_001   v1.0 simpleParameter		
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: Lack of fluidity in movement, Fluid movement,		

#### Skin texture - snout TCPIP\_CSD\_063\_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true Options: As expected, Not as expected, Locomotor activity TCPIP\_CSD\_032\_001 | v1.2 simpleParameter Req. Analysis: false Req. Upload: true Is Annotated: true Unit Measured: Squares crossed Hindlimbs - size TCPIP\_CSD\_022\_001 | v1.1 simpleParameter Req. Analysis: false Req. Upload: true Is Annotated: true Options: As expected, Both hindlimbs - decreased, One hindlimb only - decreased, Both hindlimbs - increased, One hindlimb only - increased,

Coat - color pattern - abdomen TCPIP\_CSD\_010\_001 | v1.0

Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: As expected, No dat	a, Not as expected,	
Forelimb digit - nul	mber TCPIP_CSD_044_00	2   v2.0
Req. Analysis: false	Req. Upload: true	Is Annotated: true
•	ased number, One forelimb - dember, As expected, Both foreling	
Size of squares in a procedureMetadata	arena TCPIP_CSD_086_0	01   v1.0
Req. Analysis: true	Req. Upload: true	Is Annotated: false
Unit Measured: cm^2		

# Forelimb digit - shape TCPIP\_CSD\_048\_002 | v2.0

simpleParameter

Options: One forelimb - not as expected, Both forelimbs - not as expected, As expected,		
Activity (body posi simpleParameter	tion) TCPIP_CSD_029_00	01   v1.0
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: Increased activity, In	active, As expected,	
Head morphology simpleParameter	TCPIP_CSD_026_001   v1.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: Not as expected, As	expected,	
Transfer arousal TCPIP_CSD_031_001   v1.2 simpleParameter		
Req. Analysis: false	Req. Upload: false	Is Annotated: false
Options: As expected, Immediate movement, Extended freeze,		

# Coat - color - back TCPIP\_CSD\_005\_001 | v1.0

simpleParameter

Reg. Analysis: false Reg. Upload: true Is Annotated: true **Options:** As expected, No data, Not as expected, Coat - hair distribution - abdomen TCPIP\_CSD\_013\_001 | v1.0 simpleParameter Req. Analysis: false Req. Upload: true Is Annotated: true Options: As expected, Not as expected, Genitalia - morphology TCPIP\_CSD\_073\_001 | v1.1 simpleParameter Req. Analysis: false Req. Upload: true Is Annotated: true Options: Not as expected, No data, As expected,

Vibrissae - appearance TCPIP\_CSD\_024\_001 | v1.0

Req. Analysis: false Req. Upload: true Is Annotated: true **Options:** Not as expected, No data, As expected, Hindlimb nail - shape TCPIP\_CSD\_055\_002 | v2.0 simpleParameter Reg. Analysis: false Reg. Upload: true Is Annotated: true Options: No data, Both hindlimbs - not as expected, As expected, One hindlimb - not as expected, Forelimb nail - length TCPIP\_CSD\_052\_002 | v2.0 simpleParameter Req. Analysis: false Req. Upload: true Is Annotated: true Options: No data, One forelimb - decreased length, Both forelimbs - increased length, One forelimb - increased length, Both forelimbs - decreased length, As expected,

#### Head bobbing TCPIP\_CSD\_080\_001 | v1.0

simpleParameter

Options: Present, Absent,		
Skin texture - ear TesimpleParameter	CPIP_CSD_064_001   v1.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: Not as expected, As	expected,	
Coat - hair distribution - back TCPIP_CSD_011_001   v1.0 simpleParameter		
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: As expected, Not as	expected,	
Coat - color pattern - back TCPIP_CSD_008_001   v1.0 simpleParameter		
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: As expected, No dat	ta, Not as expected,	

#### Hindlimb digit - size TCPIP\_CSD\_047\_001 | v1.2

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

**Options:** One hindlimb - increased size, As expected, Both hindlimbs - increased size, Both hindlimbs - decreased size, One hindlimb - fused, Both hindlimbs - fused,

One hindlimb - decreased size,

#### Head size TCPIP\_CSD\_025\_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

**Options:** Increased, Decreased, As expected,

# Skin color - whole body TCPIP\_CSD\_056\_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: Not as expected, As expected,

.....

Skin texture - front paws TCPIP\_CSD\_065\_001 | v1.0

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: Not as expected, As expected,

.....

#### Trunk curl TCPIP\_CSD\_038\_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: No data, Present, Absent,

.....

# Tail - morphology TCPIP\_CSD\_004\_001 | v1.1

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

**Options:** As expected, No data, Not as expected,

#### Snout size TCPIP\_CSD\_028\_001 | v1.0

simpleParameter

Options: Decreased, As expected, Increased,		
Forelimb nail - sha	<b>PE</b> TCPIP_CSD_054_002	v2.0
Req. Analysis: false	Req. Upload: true	Is Annotated: true
<b>Options:</b> No data, Both forelin As expected,	nbs - not as expected, One fore	elimb - not as expected,
Startle response TO simpleParameter	CPIP_CSD_036_001   v1.1	
Req. Analysis: false	Req. Upload: false	Is Annotated: true
Options: None, Present, No d	ata,	
Forepaw - shape TO simpleParameter	CPIP_CSD_042_002   v2.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: One forepaw - not as	s expected, As expected, Both	forepaws - not as expected,

\_\_\_\_\_\_

Tail elevation TCPIP_CSD_034_001   v1.2 simpleParameter			
Req. Analysis: false	Req. Upload: false	Is Annotated: false	
Options: As expected, No date	a, Straub / elevated tail, Dragg	ing,	
Coat - hair distribution - head TCPIP_CSD_012_001   v1.0 simpleParameter			
Req. Analysis: false	Req. Upload: true	Is Annotated: true	
Options: As expected, Not as expected,			
Hindlimb nail - number TCPIP_CSD_051_002   v2.0 simpleParameter			
Req. Analysis: false	Req. Upload: true	Is Annotated: true	
Options: Both hindlimbs - decreased number, As expected, Both hindlimbs - increased number, One hindlimb - decreased number, No data, One hindlimb - increased number,			

------

#### Upper lip morphology TCPIP\_CSD\_075\_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true **Options:** As expected, Not as expected, Coat - hair texture / appearance - abdomen TCPIP\_CSD\_016\_001 I v1.0 simpleParameter Req. Analysis: false Req. Upload: true Is Annotated: true Options: Not as expected, As expected, No data, Hindlimbs - shape TCPIP\_CSD\_021\_001 | v1.1 simpleParameter Req. Analysis: false Req. Upload: true Is Annotated: true Options: As expected, Both hindlimbs - not as expected, One hindlimb only - not as expected,

Forelimbs - shape TCPIP\_CSD\_018\_001 | v1.0

Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: One forelimb only -	not as expected, Both forelimbs	s - not as expected, As expected,
Skin color - front p simpleParameter	aws TCPIP_CSD_059_001	I   v1.0
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: As expected, Not as	expected,	
Experimenter ID TO	PIP_CSD_081_001   v1.0	
procedureMetadata		
Req. Analysis: false	Req. Upload: true	Is Annotated: false
Forelimbs - position TCPIP_CSD_017_001   v1.0		
simpleParameter		
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: As expected, One for	orelimb only - not as expected,	Both forelimbs - not as expected,

# Genitalia - presence TCPIP\_CSD\_071\_001 | v1.1

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true Options: As expected, Absent, Vibrissae - presence TCPIP\_CSD\_023\_001 | v1.0 simpleParameter Req. Analysis: false Req. Upload: true Is Annotated: true Options: As expected, Not as expected, Location of test TCPIP\_CSD\_082\_001 | v1.1 procedureMetadata Req. Analysis: false Req. Upload: true Is Annotated: false Options: LAF cabinet, Open bench,

Vocalization TCPIP\_CSD\_078\_001 | v1.0

Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: Not as expected, As	s expected,	
Number of days sin	nce cage changed 1	TCPIP_CSD_084_001   v1.1
Req. Analysis: false	Req. Upload: false	Is Annotated: false
Mouth morphology simpleParameter	TCPIP_CSD_074_001   v1.	.0
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: Not as expected, As	s expected,	
Upper teeth appearance TCPIP_CSD_069_001   v1.0 simpleParameter		
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: As expected, No date	ta, Not as expected,	

\_\_\_\_\_\_

# Comments (in English) TCPIP\_CSD\_087\_001 | v1.2

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

#### Skin color - snout TCPIP\_CSD\_057\_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: As expected, Not as expected,

#### Tail - thickness TCPIP\_CSD\_003\_001 | v1.1

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: No data, Decreased, Increased, As expected,

.....

# Hindlimb nail - length TCPIP\_CSD\_053\_002 | v2.0

simpleParameter

Options: As expected, No data, One hindlimb - increased length,

One hindlimb - decreased length, Both hindlimbs - increased length,

Both hindlimbs - decreased length,

#### Tail - length TCPIP\_CSD\_002\_001 | v1.1

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: Short, As expected, No data, Long,

# Forelimb digit - size TCPIP\_CSD\_046\_001 | v1.1

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: One forelimb - decreased size, Both forelimbs - fused,

One forelimb - increased size, One forelimb - fused, Both forelimbs - increased size,

As expected, Both forelimbs - decreased size,

#### Contact righting TCPIP\_CSD\_077\_001 | v1.0

simpleParameter

Options: Present, Absent, No data,		
Ears TCPIP_CSD_027_00 simpleParameter	01   v1.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: As expected, Not as	expected,	
Skin color - ear TCP simpleParameter	IP_CSD_058_001   v1.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: As expected, Not as	expected,	
Lower lip morpholosimpleParameter	Dgy TCPIP_CSD_076_001	v1.0
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: As expected, Not as	expected,	

# Coat - hair texture / appearance - back TCPIP\_CSD\_014\_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: No data, As expected, Not as expected,

#### Genitalia - size TCPIP\_CSD\_072\_001 | v1.1

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: As expected, No data, Not as expected,

.....

# Hindpaw - size TCPIP\_CSD\_041\_002 | v2.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

**Options:** Both hindpaws - decreased, Both hindpaws - increased, One hindpaw - increased, As expected, One hindpaw - decreased,

.....

Coat - color - abdomen TCPIP\_CSD\_007\_001 | v1.0

**Req. Analysis:** false **Req. Upload:** true **Is Annotated:** true **Options:** No data, Not as expected, As expected, Hindlimb digit - shape TCPIP\_CSD\_049\_002 | v2.0 simpleParameter Req. Analysis: false Req. Upload: true Is Annotated: true Options: One hindlimb - not as expected, As expected, Both hindlimbs - not as expected, Images TCPIP\_CSD\_085\_001 | v1.1 seriesMediaParameter Req. Analysis: false Req. Upload: false Is Annotated: false Hindlimbs - position TCPIP\_CSD\_020\_001 | v1.1 simpleParameter Req. Analysis: false Req. Upload: true Is Annotated: true

**Options:** One hindlimb only - not as expected, Both hindlimbs - not as expected, As expected,

Skin texture - whole body TCPIP_CSD_062_001   v1.0 simpleParameter		
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: Not as expected, As	expected,	
Forepaw - size TCPIF simpleParameter	P_CSD_040_002   v2.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: As expected, One forepaw - decreased, One forepaw - increased, Both forepaws - decreased, Both forepaws - increased,		
Hindpaw - shape TCPIP_CSD_043_002   v2.0 simpleParameter		
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: One hindpaw - not as expected, Both hindpaws - not as expected, As expected,		

#### Lower teeth appearance TCPIP\_CSD\_070\_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: As expected, No data, Not as expected,

#### Coat - color pattern - head TCPIP\_CSD\_009\_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: As expected, Not as expected, No data,

# Number of animals in cage TCPIP\_CSD\_083\_001 | v1.2

procedureMetadata

Req. Analysis: false Req. Upload: true Is Annotated: false

.....

# Skin texture - tail TCPIP\_CSD\_067\_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

**Options:** As expected, Not as expected, No data,

Limb grasp TCPIP_CS	SD_039_001   v1.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: No data, Present, Ab	osent,	
Skin color - back pasimpleParameter	<b>aws</b> TCPIP_CSD_060_001	1   v1.0
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: Not as expected, As	expected,	
Aggression TCPIP_CS	SD_079_001   v1.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: Excessive aggressio	n, As expected,	

#### Teeth presence TCPIP\_CSD\_068\_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: As expected, Not as expected,

#### Forelimbs - size TCPIP\_CSD\_019\_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: Both forelimbs - increased, As expected, One forelimb only - increased,

One forelimb only - decreased, Both forelimbs - decreased,

.....

#### Tail - presence TCPIP\_CSD\_001\_001 | v1.1

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Options: Absent, Present,

.....

#### Forelimb nail - number TCPIP CSD 050 002 | v2.0

simpleParameter

Req. Analysis: false	Req. Upload: true	Is Annotated: true
	orelimbs - decreased number, Noer, Both forelimbs - increased liber,	
Hindlimb digit - nu	mber TCPIP_CSD_045_00	02   v2.0
Req. Analysis: false	Req. Upload: true	Is Annotated: true
•	ased number, Both hindlimbs - ecreased number, Both hindlim	
Skin color - tail TCP simpleParameter	IP_CSD_061_001   v1.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: Not as expected, As	expected, No data,	

Skin texture - back paws TCPIP\_CSD\_066\_001 | v1.0

simpleParameter

Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: Not as expected, As expected,		
Tremor TCPIP_CSD_03 simpleParameter	0_001   v1.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Options: Present, Absent,		

\_\_\_\_\_\_