

Indirect ophthalmoscopy HAS_OPH_001

Parameters and Metadata

Eye HAS_OPH_001_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: Eye

Options: normal, abnormal,

Description HAS_OPH_002_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: Description

Left retina HAS_OPH_003_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Retina

Options: normal, abnormal, no data,

Left retina pigmentation HAS_OPH_004_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Retina_Pigmentation

Options: normal, abnormal, not defined, no data,

Left retina structure HAS_OPH_005_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Retina_Structure

Options: normal, degenerate, not defined, no data,

Left optic disc HAS_OPH_006_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Optic_Disc

Options: normal, abnormal, not defined, no data,

Right retina HAS_OPH_007_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Retina

Options: normal, abnormal, no data,

Right retina pigmentation HAS_OPH_008_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Retina_Pigmentation

Options: normal, abnormal, not defined, no data,

Right retina structure HAS_OPH_009_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Retina_Structure

Options: normal, degenerate, not defined, no data,

Right optic disc HAS_OPH_010_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Optic_Disc

Options: normal, abnormal, not defined, no data,

Left eye blood vessels HAS_OPH_011_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Eye_Blood_Vessels

Options: normal, abnormal, no data,

Left eye blood vessels structure HAS_OPH_012_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Eye_Blood_Vessels_Structure

Options: normal appearance, abnormal appearance, not defined, no data,

Left eye blood vessels pattern HAS_OPH_013_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Eye_Blood_Vessels_Pattern

Options: normal, abnormal, not defined, no data,

Right eye blood vessels HAS_OPH_014_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Eye_Blood_Vessels

Options: normal, abnormal, no data,

Right eye blood vessels structure HAS_OPH_015_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Eye_Blood_Vessels_Structure

Options: normal appearance, abnormal appearance, not defined, no data,

Right eye blood vessels pattern HAS_OPH_016_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Eye_Blood_Vessels_Pattern

Options: normal, abnormal, not defined, no data,

Digital image(s) HAS_OPH_017_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: Digital_Image_s_

Increments: Minimum 1

Equipment name HAS_OPH_018_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_name

Equipment manufacturer HAS_OPH_019_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_manufacturer

Equipment model HAS_OPH_020_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_model

Body weight HAS_OPH_021_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

Description: Body_Weight

Topical agent(s) HAS_OPH_022_001 | v1.0

procedureMetadata

Req. Analysis: true **Req. Upload:** true **Is Annotated:** false

Description: Topical_agent_s_

Options: atropine, tropicamide, atropine + neosynephrine, tropicamide + neosynephrine, tropicamide + phenylephrine,