

# X-Ray HAROLD\_006\_001

## Parameters and Metadata

### Clavicle HAROLD\_006\_001\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: Clavicle

Options: imageOnly, normal, abnormal,

---

### Completeness of digits HAROLD\_006\_001\_002 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: Completeness\_of\_digits

Options: imageOnly, yes, no,

---

### Digit integrity HAROLD\_006\_001\_003 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

**Description:** Digit\_integrity

**Options:** imageOnly, yes, no,

---

**Number of digits** HAROLD\_006\_001\_004 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Number\_of\_digits

---

**Shape of femur** HAROLD\_006\_001\_005 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Shape\_of\_femur

**Options:** imageOnly, normal, abnormal,

---

**Fibula** HAROLD\_006\_001\_006 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Fibula

**Options:** imageOnly, normal, abnormal,

---

**Shape of humerus** HAROLD\_006\_001\_007 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Shape\_of\_humerus

**Options:** imageOnly, normal, abnormal,

---

**Joints** HAROLD\_006\_001\_008 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Joints

**Options:** imageOnly, normal, abnormal,

---

**Mandibles** HAROLD\_006\_001\_009 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Mandibles

**Options:** imageOnly, normal, abnormal,

---

**Maxilla** HAROLD\_006\_001\_010 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Maxilla

**Options:** imageOnly, normal, abnormal,

---

**Orbit** HAROLD\_006\_001\_011 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Orbit

**Options:** imageOnly, normal, abnormal,

---

**Pelvis** HAROLD\_006\_001\_012 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Pelvis

**Options:** imageOnly, normal, abnormal,

---

**Radius** HAROLD\_006\_001\_013 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Radius

**Options:** imageOnly, normal, abnormal,

---

**Rib** HAROLD\_006\_001\_014 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Rib

**Options:** imageOnly, fused, not fused,

---

**Number of ribs left** HAROLD\_006\_001\_015 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: Number\_of\_ribs\_left

---

**Number of ribs right** HAROLD\_006\_001\_016 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: Number\_of\_ribs\_right

---

**Shape of ribs** HAROLD\_006\_001\_017 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: Shape\_of\_ribs

Options: imageOnly, normal, abnormal,

---

**Scapulae** HAROLD\_006\_001\_018 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

**Description:** Scapulae

**Options:** imageOnly, normal, abnormal,

---

**Skull shape** HAROLD\_006\_001\_019 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Skull\_shape

**Options:** imageOnly, normal, abnormal,

---

**Syndactylism** HAROLD\_006\_001\_020 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Syndactylism

**Options:** imageOnly, syndactylism, no syndactylism,

---

**Teeth** HAROLD\_006\_001\_021 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Teeth

**Options:** imageOnly, normal, abnormal,

---

**Tibia** HAROLD\_006\_001\_022 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Tibia

**Options:** imageOnly, normal, abnormal,

---

**Ulna** HAROLD\_006\_001\_023 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Ulna

**Options:** imageOnly, normal, abnormal,

---

**Number of caudal vertebrae** HAROLD\_006\_001\_024 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true



**Description:** Number\_of\_caudal\_vertebrae

---

**Number of cervical vertebrae** HAROLD\_006\_001\_025 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** true

**Description:** Number\_of\_cervical\_vertebrae

---

**Fusion of vertebrae** HAROLD\_006\_001\_026 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** true

**Description:** Fusion\_of\_vertebrae

**Options:** imageOnly, yes, no,

---

**Number of lumbar vertebrae** HAROLD\_006\_001\_027 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** true

**Description:** Number\_of\_lumbar\_vertebrae

---

## Number of pelvic vertebrae HAROLD\_006\_001\_028 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Number\_of\_pelvic\_vertebrae

---

## Processes on vertebrae HAROLD\_006\_001\_029 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Processes\_on\_vertebrae

**Options:** imageOnly, normal, abnormal,

---

## Shape of vertebrae HAROLD\_006\_001\_030 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Shape\_of\_vertebrae

**Options:** imageOnly, normal, abnormal,

---

## Number of thoracic vertebrae HAROLD\_006\_001\_031 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Description:** Number\_of\_thoracic\_vertebrae

---

## Comment on X-Ray descriptive part HAROLD\_006\_001\_032 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

**Description:** Comment\_on\_XRay\_descriptive\_part

---

## X-Ray images HAROLD\_006\_001\_033 | v1.0

seriesMediaParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Description:** XRay\_images

**Increments:** Minimum 1

---

## **Alive** HAROLD\_006\_001\_034 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Description:** Alive\_or\_Dead

**Options:** yes, no,

---

## **Tibia Length** HAROLD\_006\_001\_035 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Unit Measured:** mm

**Description:** Tibia\_Length

---

## **Radio-opacities** HAROLD\_006\_001\_036 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

**Description:** radioopacities

**Options:** No, Yes,

---

## Date of procedure HAROLD\_006\_001\_037 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: date\_of\_procedure

---

## Xray Equipment name HAROLD\_006\_001\_801 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Xray\_Equipment\_name

---

## Xray Equipment manufacturer HAROLD\_006\_001\_802 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Xray\_Equipment\_manufacturer

---

## Xray Equipment model HAROLD\_006\_001\_803 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Xray\_Equipment\_model

---

## Settings time of exposure HAROLD\_006\_001\_804 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Settings\_time\_of\_exposure

---

## Settings level HAROLD\_006\_001\_805 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Settings\_level

---

## Experimenter (Scan) HAROLD\_006\_001\_806 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: Experimenter\_Scan

---

## Experimenter (Analysis) HAROLD\_006\_001\_807 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Experimenter\_Analysis

---

## Voltage settings HAROLD\_006\_001\_808 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Voltage\_settings

---

## Scanner Equipment name HAROLD\_006\_001\_809 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: Scanner\_Equipment\_name

---

## Scanner Equipment manufacturer HAROLD\_006\_001\_810 | v1.0

procedureMetadata

Req. Analysis: false      Req. Upload: false      Is Annotated: false

Description: Scanner\_Equipment\_manufacturer

---

**Scanner Equipment model** HAROLD\_006\_001\_811 | v1.0

procedureMetadata

Req. Analysis: false      Req. Upload: false      Is Annotated: false

Description: Scanner\_Equipment\_model

---