

Clinical Chemistry HAROLD_015_001

Parameters and Metadata

Glucose HAROLD_015_001_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mmol/l

Description: Glucose

Urea HAROLD_015_001_002 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mmol/l

Description: Urea

Creatinine HAROLD_015_001_003 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: umol/l

Description: Creatinine

Sodium HAROLD_015_001_004 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mmol/l

Description: Sodium

Potassium HAROLD_015_001_005 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mmol/l

Description: Potassium

Chloride HAROLD_015_001_006 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mmol/l

Description: Chloride

Total protein HAROLD_015_001_007 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: g/l

Description: Total_protein

Albumin HAROLD_015_001_008 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: g/l

Description: Albumin

Calcium HAROLD_015_001_009 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mmol/l

Description: Calcium

Phosphorus HAROLD_015_001_010 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mmol/l

Description: Phosphorus

Iron HAROLD_015_001_011 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: umol/l

Description: Iron

Lactate dehydrogenase HAROLD_015_001_012 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: U/l

Description: Lactate_dehydrogenase

Aspartate aminotransferase HAROLD_015_001_013 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: U/l

Description: Aspartate_aminotransferase

Alanine aminotransferase HAROLD_015_001_014 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: U/l

Description: Alanine_aminotransferase

Alkaline phosphatase HAROLD_015_001_015 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: U/l

Description: Alkaline_phosphatase

Alpha-amylase HAROLD_015_001_016 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: U/l

Description: Alphaamylase

Total cholesterol HAROLD_015_001_017 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mmol/l

Description: Total_cholesterol

Triglyceride HAROLD_015_001_018 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mmol/l

Description: Triglyceride

Free fatty acid HAROLD_015_001_019 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mmol/l

Description: Free_fatty_acid

Creatine kinase HAROLD_015_001_020 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: U/l

Description: Creatine_kinase

Uric acid HAROLD_015_001_021 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: umol/l

Description: Uric_acid

Total bilirubin HAROLD_015_001_022 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: umol/l

Description: Total_bilirubin

HDL-cholesterol HAROLD_015_001_023 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mmol/l

Description: HDLcholesterol

LDL-cholesterol HAROLD_015_001_024 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mmol/l

Description: LDLcholesterol

Ferretin HAROLD_015_001_025 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mg/ml

Description: Ferretin

Transferrin HAROLD_015_001_026 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mg/dl

Description: Transferrin

C-reactive protein HAROLD_015_001_027 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mg/l

Description: CReactive_protein

Fructosamine HAROLD_015_001_028 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: umol/l

Description: fructosamine

Equipment name HAROLD_015_001_801 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: Equipment_name

Equipment manufacturer HAROLD_015_001_802 | v1.0

procedureMetadata

Req. Analysis: true Req. Upload: true Is Annotated: false

Description: Equipment_manufacturer

Equipment model HAROLD_015_001_803 | v1.0

procedureMetadata

Req. Analysis: true Req. Upload: true Is Annotated: false

Description: Equipment_model

Method of blood collection HAROLD_015_001_804 | v1.0

procedureMetadata

Req. Analysis: true Req. Upload: true Is Annotated: false

Description: Method_of_blood_collection

EMPreSSID for blood collection SOP HAROLD_015_001_805 | v1.0

procedureMetadata

Req. Analysis: false Req. Upload: false Is Annotated: false

Description: EMPReSSID_for_blood_collection_SOP

Date/time of blood collection HAROLD_015_001_806 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: DateTime_of_blood_collection

Fasting prior to experiment HAROLD_015_001_807 | v1.0

procedureMetadata

Req. Analysis: true **Req. Upload:** false **Is Annotated:** false

Description: Fasting_prior_to_experiment

Options: yes, no,

Approximate period of fasting HAROLD_015_001_808 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Unit Measured: Hours

Description: Period_of_fasting

Moved from cage for fasting HAROLD_015_001_809 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Description: Moved_from_cage_for_fasting

Options: yes, no,

Plasma dilution HAROLD_015_001_810 | v1.0

procedureMetadata

Req. Analysis: true **Req. Upload:** true **Is Annotated:** false

Description: Plasma_dilution

Sample status HAROLD_015_001_811 | v1.0

procedureMetadata

Req. Analysis: true **Req. Upload:** true **Is Annotated:** false

Description: Sample_Status

Options: fresh, frozen,

Anaesthesia used for blood collection HAROLD_015_001_812 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: Anaesthesia_used_for_blood_collection

Date of measurement HAROLD_015_001_813 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Day_of_measurement

Haemolysis status HAROLD_015_001_814 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: haemolysis_status

Options: N, /, NP, +, ++, +++, +++++, ++++++, ABN,
