

Body Weight RBRCIP_BWT_001

Purpose

The body weight test measures the weight of the mouse in a time series, allowing monitoring of its evolution; also, it is required in many other procedures.

Equipment

- Laboratory Balance

Procedure

- Place a cage containing mice on an operation table
- Remove a mouse from its cage. Determine and record its weight.
- Return the mouse to the cage or continue with a procedure.

Notes

The body weight of each mouse is to be measured every two months, from the 24th until the 48th week. Optional additional weights may also be recorded.

Parameters and Metadata

Equipment manufacturer RBRCIP_BWT_004_001 | v1.2

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Options: Scientech, Denver Instrument, Acculab, Ohaus, Mettler Toledo, Kern & Sohn GmbH, Sartorius, Radwag,

Body weight curve RBRCIP_BWT_008_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Increments: Minimum 1

Equipment model RBRCIP_BWT_007_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Options: CPA3202S, PCB2000-1, EMB220-1, BL310, SL-3100D, EMB500-1, TP-202, BP6100, VIC-511, 440-47, EMB 200-2, AM100, EMB600-2, 440-33, Entris 2201-1S, SI-2002, Adventurer Pro, PG3001-S, EW600-2M, 440-47N, AV213C, 770-14, AV212C, Navigator 34120, AV2101, AB104-S, VIC-123, 572-35, XS802S, 201-10,

Date equipment last calibrated RBRCIP_BWT_006_001 | v1.2

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Body weight RBRCIP_BWT_001_001 | v1.3

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

General comments about the mouse RBRCIP_BWT_002_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Equipment ID RBRCIP_BWT_003_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Experimenter ID RBRCIP_BWT_005_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false
