Sleep: 24LD activity HAS_SLE_001

Parameters and Metadata

Distance travelled simpleParameter	HAS_SLE_001_001 v1.3	
Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: cm		
Description: distance_travel	led	
Date of procedure simpleParameter	HAS_SLE_002_001 v1.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: false
Procedural comments HAS_SLE_003_001 v1.0 simpleParameter		

General comments about the mouse HAS_SLE_004_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Description: general_comments_about_the_mouse

Equipment ID HAS_SLE_005_001 | v1.1

procedureMetadata

Req. Analysis: false Req. Upload: true Is Annotated: false

Description: equipment_name

.....

Equipment manufacturer HAS_SLE_006_001 | v1.0

procedureMetadata

Req. Analysis: true Req. Upload: true Is Annotated: false

Description: equipment_manufacturer

Equipment model HAS_SLE_007_001 | v1.0

procedureMetadata

Req. Analysis: true Req. Upload: true Is Annotated: false

Coffin number HAS_SLE_008_001 | v1.2 procedureMetadata Req. Analysis: true Req. Upload: true Is Annotated: false **Options:** 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, Apparatus HAS_SLE_009_001 | v1.1 procedureMetadata Req. Analysis: true Req. Upload: true Is Annotated: false Mean speed HAS_SLE_011_001 | v1.2 simpleParameter Req. Analysis: false Req. Upload: false Is Annotated: true

Mobile episodes HAS_SLE_014_001 | v1.3

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Immobile episodes HAS_SLE_015_001 | v1.3

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Time segment of test HAS_SLE_016_001 | v1.2

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Unit Measured: s

Options: 0 - 3600, 3600 - 7200, 7200 - 10800, 10800 - 14400, 14400 - 18000, 18000 - 21600, 21600 - 25200, 25200 - 28800, 28800 - 32400, 32400 - 36000,

36000 - 39600, 39600 - 43200,

Time immobile D HAS_SLE_012_002 | v2.2

seriesParameter

Req. Analysis: false Req. Upload: true Is Annotated: false

Unit Measured: s

Increments: 0, 3600, 7200, 10800, 14400, 18000, 21600, 25200, 28800, 32400, 36000, 39600,

Time immobile L HAS_SLE_013_002 | v2.2

seriesParameter

Req. Analysis: false Req. Upload: true Is Annotated: false

Unit Measured: s

Increments: 0, 3600, 7200, 10800, 14400, 18000, 21600, 25200, 28800, 32400, 36000,

39600,

Time immobile HAS_SLE_010_001 | v1.1

seriesParameter

Req. Analysis: true Req. Upload: false Is Annotated: false

Unit Measured: s

Increments: 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11,

............

Immobility % HAS_SLE_017_001 | v1.3

seriesParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Unit Measured: %

Description: immobility %

Increments: dark, light,
Derivation: archived('sum(immobility)/time')