

ERG (Electroretinogram) (GMC) GMC_912_001

Parameters and Metadata

A-Wave 500 Cd/m² left eye GMC_912_001_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: uV

Description: AWave_500_Cdm_left_eye

A-Wave 500 Cd/m² right eye GMC_912_001_002 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: uV

Description: AWave_500_Cdm_right_eye

A-Wave 12500 Cd/m² left eye GMC_912_001_003 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: uV

Description: AWave_12500_Cdm_left_eye

A-Wave 12500 Cd/m^Å² right eye GMC_912_001_004 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: uV

Description: AWave_12500_Cdm_right_eye

B-Wave 500 Cd/m^Å² left eye GMC_912_001_005 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: uV

Description: BWave_500_Cdm_left_eye

B-Wave 500 Cd/m^Å² right eye GMC_912_001_006 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: uV

Description: BWave_500_Cdm_right_eye

B-Wave 12500 Cd/m^Å² left eye GMC_912_001_007 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: uV

Description: BWave_12500_Cdm_left_eye

B-Wave 12500 Cd/m^Å² right eye GMC_912_001_008 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: uV

Description: BWave_12500_Cdm_right_eye

Comment on the ERG examination GMC_912_001_009 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: Comment_on_the_ERG_examination

Latency A-wave 500 Cd/m^Å² right eye GMC_912_001_010 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: ms

Description: Latency_Awave_500_Cdm_right_eye

Latency B-wave 500 Cd/m^Å² right eye GMC_912_001_011 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: ms

Description: Latency_Bwave_500_Cdm_right_eye

Latency A-wave 500 Cd/m^Å² left eye GMC_912_001_012 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: ms

Description: Latency_Awave_500_Cdm_left_eye

Latency B-wave 500 Cd/m^Å² left eye GMC_912_001_013 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: ms

Description: Latency_Bwave_500_Cdm_left_eye

Latency A-wave 12500 Cd/m^Å² left eye GMC_912_001_014 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: ms

Description: Latency_Awave_12500_Cdm_left_eye

Latency A-wave 12500 Cd/m^Å² right eye GMC_912_001_015 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: ms

Description: Latency_Awave_12500_Cdm_right_eye

Latency B-wave 12500 Cd/m² left eye GMC_912_001_016 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: ms

Description: Latency_Bwave_12500_Cdm_left_eye

Latency B-wave 12500 Cd/m² right eye GMC_912_001_017 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: ms

Description: Latency_Bwave_12500_Cdm_right_eye

Body weight at ERG examination GMC_912_001_018 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: g

Description: Body_weight_at_ERG_examination
