

Hot Plate M-G-P_012_001

Parameters and Metadata

Time of first response M-G-P_012_001_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Unit Measured: s

Description: Time_of_first_response_

Type of response M-G-P_012_001_002 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: false

Description: Type_of_response

Options: shake, lick, jump, locomotion, no reaction, other,

Comment on other behaviours M-G-P_012_001_004 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Description: Comment_on_other_behaviours

Body Weight M-G-P_012_001_005 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** true

Unit Measured: g

Description: Body_Weight

Equipment name M-G-P_012_001_801 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: Equipment_name

Equipment manufacturer M-G-P_012_001_802 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: Equipment_manufacturer

Equipment model M-G-P_012_001_803 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_model

Type of recording M-G-P_012_001_804 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Type_of_recording

Options: automated, manual,
